

第 000035 号

令和 3 年 3 月 5 日

お客様各位

ペルモビール株式会社

一部ロホクッションカバーの素材変更に関するお知らせ

拝啓

貴社ますますご清栄のこととお慶び申し上げます。平素は格別のお引立てを賜り厚くお礼申し上げます。

さて、下記の通り一部ロホクッションカバーにおいて、耐久性向上を目的とした素材の変更を実施することとなりましたのでお知らせいたします。

なお、2021 年 3 月 14 日までを段階的な変更期間とさせていただき、2021 年 3 月 15 日米国出荷分より切り替えさせていただきます。

お客様にはご不便・ご迷惑をお掛け致しますが、何卒よろしくご厚意申し上げます。今後ともロホ製品の変わらぬご愛顧の程よろしくご厚意申し上げます。

敬具

記

対象製品：ハイブリッドエリート
ヘビーデューティー

カバー素材（下層）変更前：熱可塑性ポリウレタン（TPU）を施したポリエステルニット
カバー素材（下層）変更後：オレフィン系熱可塑性エラストマー（TPO）

● ロボ・ハイブリッドエリートカバーリスト

| 品番 | 製品名 |
|------------|--------------------------|
| COV-HE1414 | COVER, HYBRID ELITE 1414 |
| COV-HE1415 | COVER, HYBRID ELITE 1415 |
| COV-HE1416 | COVER, HYBRID ELITE 1416 |
| COV-HE1515 | COVER, HYBRID ELITE 1515 |
| COV-HE1516 | COVER, HYBRID ELITE 1516 |
| COV-HE1517 | COVER, HYBRID ELITE 1517 |
| COV-HE1615 | COVER, HYBRID ELITE 1615 |
| COV-HE1616 | COVER, HYBRID ELITE 1616 |
| COV-HE1617 | COVER, HYBRID ELITE 1617 |
| COV-HE1618 | COVER, HYBRID ELITE 1618 |
| COV-HE1620 | COVER, HYBRID ELITE 1620 |
| COV-HE1716 | COVER, HYBRID ELITE 1716 |
| COV-HE1717 | COVER, HYBRID ELITE 1717 |
| COV-HE1718 | COVER, HYBRID ELITE 1718 |
| COV-HE1720 | COVER, HYBRID ELITE 1720 |
| COV-HE1816 | COVER, HYBRID ELITE 1816 |
| COV-HE1817 | COVER, HYBRID ELITE 1817 |
| COV-HE1818 | COVER, HYBRID ELITE 1818 |
| COV-HE1820 | COVER, HYBRID ELITE 1820 |
| COV-HE2017 | COVER, HYBRID ELITE 2017 |
| COV-HE2018 | COVER, HYBRID ELITE 2018 |
| COV-HE2020 | COVER, HYBRID ELITE 2020 |
| COV-HE2218 | COVER, HYBRID ELITE 2218 |
| COV-HE2220 | COVER, HYBRID ELITE 2220 |
| COV-HE2222 | COVER, HYBRID ELITE 2222 |
| COV-HE2418 | COVER, HYBRID ELITE 2418 |
| COV-HE2420 | COVER, HYBRID ELITE 2420 |

● ヘビーデューティー用カバーリスト

| 品番 | 製品名 |
|--------------|--------------------------------|
| COV-HD1010 | COVER, HEAVY DUTY 10X10 FOR MP |
| COV-HD1010LP | COVER, HEAVY DUTY 10X10 LP |
| COV-HD1011 | COVER, HEAVY DUTY 10X11 FOR MP |
| COV-HD1011LP | COVER, HEAVY DUTY 10X11 LP |
| COV-HD1012 | COVER, HEAVY DUTY 10X12 |
| COV-HD1013 | COVER, HEAVY DUTY 10X13 |
| COV-HD1014 | COVER, HEAVY DUTY 10X14 FOR MP |
| COV-HD108 | COVER, HEAVY DUTY 10X8 |
| COV-HD108LP | COVER, HEAVY DUTY 10X8 LP |
| COV-HD109 | COVER, HEAVY DUTY 10X9 FOR MP |
| COV-HD1114 | COVER, HEAVY DUTY 11X14 |
| COV-HD118 | COVER, HEAVY DUTY 11X8 |
| COV-HD119 | COVER, HEAVY DUTY 11X9 FOR MP |
| COV-HD119LP | COVER, HEAVY DUTY 11X9 LP |
| COV-HD1210 | COVER, HEAVY DUTY 12X10 FOR MP |
| COV-HD1210LP | COVER, HEAVY DUTY 12X10 LP |
| COV-HD1211 | COVER, HEAVY DUTY 12X11 FOR MP |
| COV-HD1211LP | COVER, HEAVY DUTY 12X11 LP |
| COV-HD1212 | COVER, HEAVY DUTY 12X12 FOR MP |
| COV-HD1212LP | COVER, HEAVY DUTY 12X12 LP |
| COV-HD1213 | COVER, HEAVY DUTY 12X13 |
| COV-HD1214 | COVER, HEAVY DUTY 12X14 |
| COV-HD128 | COVER, HEAVY DUTY 12X8 |
| COV-HD129 | COVER, HEAVY DUTY 12X9 |
| COV-HD1310 | COVER, HEAVY DUTY 13X10 FOR MP |
| COV-HD1310LP | COVER, HEAVY DUTY 13X10 LP |
| COV-HD1311 | COVER, HEAVY DUTY 13X11 FOR MP |
| COV-HD1311LP | COVER, HEAVY DUTY 13X11 LP |
| COV-HD1312 | COVER, HEAVY DUTY 13X12 |
| COV-HD1313 | COVER, HEAVY DUTY 13X13 |
| COV-HD1314 | COVER, HEAVY DUTY 13X14 |
| COV-HD138 | COVER, HEAVY DUTY 13X8 |
| COV-HD139 | COVER, HEAVY DUTY 13X9 |
| COV-HD1410 | COVER, HEAVY DUTY 14X10 |

| | |
|-------------|--------------------------------|
| COV-HD1411 | COVER, HEAVY DUTY 14X11 |
| COV-HD1412 | COVER, HEAVY DUTY 14X12 |
| COV-HD1413 | COVER, HEAVY DUTY 14X13 |
| COV-HD1414 | COVER, HEAVY DUTY 14X14 |
| COV-HD148 | COVER, HEAVY DUTY 14X8 |
| COV-HD149 | COVER, HEAVY DUTY 14X9 |
| COV-HD66 | HCARD COVER, HEAVY DUTY 6X6 |
| COV-HD66LP | HCARD COVER, HEAVY DUTY 6X6 L |
| COV-HD67 | HCARD COVER, HEAVY DUTY 6X7 |
| COV-HD67LP | HCARD COVER, HEAVY DUTY 6X7 L |
| COV-HD710 | COVER, HEAVY DUTY 7X10 |
| COV-HD77 | COVER, HEAVY DUTY 7X7 FOR MP,H |
| COV-HD77LP | HCARD COVER, HEAVY DUTY 7X7 L |
| COV-HD78 | COVER, HEAVY DUTY 7X8 FOR MP,H |
| COV-HD78LP | COVER, HEAVY DUTY 7X8 LP |
| COV-HD79 | COVER, HEAVY DUTY 7X9 FOR MP,H |
| COV-HD79LP | COVER, HEAVY DUTY 7X9 LP |
| COV-HD810 | COVER, HEAVY DUTY 8X10 |
| COV-HD810LP | COVER, HEAVY DUTY 8X10 LP |
| COV-HD811 | COVER, HEAVY DUTY 8X11 FOR MP |
| COV-HD87 | COVER, HEAVY DUTY 8X7 |
| COV-HD87LP | HCARD COVER, HEAVY DUTY 8X7 L |
| COV-HD88 | COVER, HEAVY DUTY 8X8 FOR MP,H |
| COV-HD88LP | COVER, HEAVY DUTY 8X8 LP |
| COV-HD89 | COVER, HEAVY DUTY 8X9 FOR MP,H |
| COV-HD89LP | COVER, HEAVY DUTY 8X9 LP |
| COV-HD910 | COVER, HEAVY DUTY 9X10 FOR MP |
| COV-HD910LP | COVER, HEAVY DUTY 9X10 LP |
| COV-HD911 | COVER, HEAVY DUTY 9X11 FOR MP |
| COV-HD911LP | COVER, HEAVY DUTY 9X11 LP |
| COV-HD912 | COVER, HEAVY DUTY 9X12 |
| COV-HD98 | COVER, HEAVY DUTY 9X8 FOR MP,H |
| COV-HD98LP | COVER, HEAVY DUTY 9X8 LP |
| COV-HD99 | COVER, HEAVY DUTY 9X9 FOR MP,H |
| COV-HD99LP | COVER, HEAVY DUTY 9X9 LP |